



Acura Digital reading micropipettes – your expectations - our tradition

來自精密工藝國度瑞士ACURA數字式微量自動吸管
重量輕巧手感極佳、具自我校正功能、可任意調整去Tip間隙、幾乎適用所有廠牌Tips、全支可高溫高壓滅菌



型號	範圍ul	讀取量ul	精度CV%	定價
825.0002	0.1 - 2	0.002	<1.05%	7000
825.0010	0.5 - 10	0.01	<0.5%	7000
825.0020	2 - 20	0.02	<0.5%	7000
825.0050	5 - 50	0.1	<0.5%	7000
825.0100	10 - 100	0.1	<0.2%	7000
825.0200	20 - 200	0.2	<0.2%	7000
825.1000	100 - 1000	1	<0.2%	7000
<hr/>				
835.05	0.5 - 5 (ml)	0.01 (ml)	<0.2%	8000
835.10	1 - 10 (ml)	0.01 (ml)	<0.2%	8000
<hr/>				
855.08.010	0.5 - 10	0.01	<1.0%	17800
855.08.050	5 - 50	0.1	<0.4%	17800
855.08.350	40 - 350	0.01	<0.3%	17800
<hr/>				
855.12.010	0.5 - 10	0.01	<1.0%	18800
855.12.050	5 - 50	0.1	<0.4%	18800
855.12.350	40 - 350	0.01	<0.3%	18800



可任意調整
去Tip間隙

SOCOREX 於1960由2位瑞士錶匠所創立

台灣轉錄科技公司

TEL: (07)5542589 FAX: (07)5553545 Toll Free: 0800 231 233

